

Using Penguin^{RFA} together with an Osstell[®] Smartpeg[®]

Measured ISQ	Add	Measured ISQ	Add	Measured ISQ	Add
20	4	44	0	68	0
21	4	45	1	69	0
22	3	46	1	70	0
23	3	47	2	71	0
24	3	48	2	72	-1
25	3	49	2	73	-1
26	3	50	2	74	-1
27	3	51	2	75	-1
28	2	52	2	76	-2
29	2	53	3	77	-2
30	2	54	3	78	-2
31	2	55	3	79	-3
32	2	56	3	80	-3
33	2	57	2	81	-3
34	2	58	2	82	-3
35	1	59	2	83	-4
36	1	60	2	84	-4
37	1	61	2	85	-4
38	1	62	2	86	-4
39	1	63	2	87	-5
40	1	64	1	88	-5
41	1	65	1	89	-5
42	1	66	1	90	-6
43	1	67	1		

Take the value in the "Add" column and add to the measured ISQ to achieve the correct ISQ value

