






Taking Your KaVo NOMAD™ Pro 2 Battery Out of Sleep Mode



Sometimes when your handset has not been used for an extended amount of time, it may enter sleep mode. Luckily, waking up the battery is quick, easy and safe.

How do you know if your battery is in sleep mode? If the handset is properly placed on the charging cradle and the lights do not illuminate, the battery may be in sleep mode. If so, simply follow the steps below to wake up the battery.

Step One	Step Two	Step Three
Plug in the charger and place the handset on the charger.	Leaving the handset on the charger, unplug the charger by unplugging the power supply from either the cradle or the wall.	Wait 10 seconds and then plug the power supply back into the cradle or the wall.
		

This will reset the handset and it will begin charging (lights on charging cradle will begin cycling through).

Note: Charging a handset from the sleep mode may take slightly longer than normal. If the process described above fails to recover the handset, authorized service may be required.

Kavo™ and NOMAD™ are trademarks of Kavo Kerr affiliates in the United States and/or other countries. The products, equipment and services illustrated and described in this piece reflect knowledge at the time of printing. KaVo Dental accepts no liability for any deviation from the illustrations in terms of color or shape, or any errors or print errors, and retains the right to make changes to the brochure at any time.

For indications for use, please visit: www.kavo.com/us/ifu.

Dental Imaging Technologies Corporation | 1910 North Penn Rd | Hatfield, PA 19440 | USA | 1-888-ASK-KAVO

www.kavo.com

© 2018 KaVo Dental. All rights reserved. KV00570/RevA