



A Portable Dental System That
SHATTERS BARRIERS
of size, weight, performance and
energy efficiency.

Introducing Aseptico's most advanced portable dental system ever – the GO Ultra Portable.

- **SUPER-COMPACT** – So small it qualifies as carry-on luggage
- **LIGHTWEIGHT** – 10 lbs lighter than the world's #1 system
- **PERFORMANCE** – New vacuum blower technology generates powerful suction and air flow
- **ULTRA-ENERGY EFFICIENT** – provides new possibilities for operating on battery power

It's a joy to take dentistry or hygiene care anywhere, with the new Aseptico GO Ultra Portable Dental System. Ideal for urban, rural, and remote dental care, the GO is capable of restorative dentistry and hygiene functions with an integrated electric motor and LED fiber-optic lighting. The GO is fully regulatory compliant, features universal voltage inputs, and is certified for electrical safety. The new GO is neatly contained in an impact resistant, waterproof wheeled case with a soft handle that is easy on the hands. It's everything a professional needs to take dentistry on the GO!



Faster Drying for Sealants
Syringe Performance:
80 psi initial
50 psi continuous

Quiet & Powerful Electric Motor
5k to 200k rpm range is ideal for prophylaxis and restorative.

More Powerful Suction!
New vacuum blower technology generates 5 SCFM

1 Liter Waste Capacity
Easy, clean removal of waste

500 mL Water Capacity
Plenty of water capacity for restorative cases

Impact Resistant

Durable, wheeled case with comfortable handles.

Super Compact!

22" x 13" x 9.2" fits into the airline overhead as carry on luggage

Easy and Light

Only 35 pounds!
10 pounds lighter than the Express, the world's #1 system!

Plug in Anywhere

100-240 VAC & 24VDC input options



Weighted and Omnidirectional Foot Pedal
Stays in place and can be activated from any angle

