

Aseptico

2025

DENTAL CATALOG

CELEBRATING

50

YEARS IN
DENTISTRY

PORTABLE DENTISTRY

NEW Scaler Option on the
GO Ultra Portable Dental System (pg. 32)

INNOVATIONS & UTILITY

Assess Osseointegration: Penguin II (pg. 22)


DENTAL MOTORS & HANDPIECES

NEW Implant Motor: 1070 Series (pg. 4)



NEW CHAIRS & STOOLS UPDATE
SILEATHER FABRIC
Sustainable & Stronger





Welcome to the 2025 Aseptico Global Dental Catalog

This year marks Aseptico’s 50th year serving the dental profession. It is a privilege to create equipment for dental professionals worldwide, whose dedication to patient care continually inspires us.

Each year, we improve our products to help clinicians perform dentistry more effectively – whether placing implants in a high-tech practice or delivering dental care in non-traditional locations such as nursing homes, school sealant programs, or humanitarian clinics.

In 2024, Aseptico introduced six new products including two new oral surgery handpieces, a new implant hand-piece, the ATC-16 Deluxe Operatory Cart, the new GO Ultra Portable AEU-350S portable dental system with built-in ultrasonic scaler, and the new GO Cart Portable Supply Cart.

Additionally, the 1070 Series Implant Motor, our flagship 7000 Series All-in-one platform motor, the AMC-15 Mobile Dental Cart, and the new GO Ultra Portable AEU-350S received significant recognition for their performance capabilities from Dentistry Today, the nation’s leading clinical news magazine for dentists.

You will find this year’s catalog pages are filled with all of these new products, plus our full line of motor systems and handpieces, innovative and trusted products such as HandiDam pre-framed dental dams, Penguin II Implant Stability Monitor, Spotter Smart Implant Detector, and more.

We take pride that almost all Aseptico products are made in the USA and invite you to compare Aseptico products to any other for quality, performance, innovative features, reliability, and value.

We look forward to serving you in the coming year and beyond!

Glenn Kazen
President of Aseptico Inc.

Aseptico DENTAL CATALOG



Dental Motors & Handpieces

WE POWER DENTISTRY

Recommended Systems	3
Oral Surgery and Implant Motors	3-8
1070 Series Implant Motor	4-5
Electrical Operative Motor	9
Endodontic Motors	10-11
NEW Handpieces	12-18
Tubing, Accessories & Replacement Parts	19-20



Innovations & Utility

PROVEN PRODUCTS

Penguin II Implant Stability Monitor	22
Spotter Smart Implant Detector	23
Powered Trolley Carts	26
NEW Deluxe Operatory Cart	26
HandiDam® Pre-Framed Dental Dam	27
Ortho Technique Chain & Grid	27
Dental Dam Clamps	27
Dental Dam Punches and Forceps	27
Dental Delivery Systems	28



Mobile & Portable Dentistry

MOVING TO MEET NEED

Recommended Systems	31
Portable Delivery Systems	31-36
GO Ultra Portable Delivery System	
(NEW Scaler Option Available)	32-33
NEW GO Cart Portable Supply Cart	34
Mobile Dental Carts	38-40
Portable Chairs & Stools	42-43
NOMAD Pro Portable X-ray	44
Portable LED Exam Light	45
Portable Tray Stand	45
Lab Motors	46
Portable Saliva Ejector	47

Customer Service & Ordering

VALUE & EXCELLENCE

Ordering Information, Policies and Contact Information	48
Repair and Product Support	48

Call 425.408.9347 or Shop Online at www.aseptico.com

Aseptico
DENTAL MOTORS
AND HANDPIECES



Find The Right Dental Motor For You

Driven to advancing dental care and patient access to care. **Call us today! 425-408-9347.**

				
	Our affordable new implant motor streamlines procedures	Award-winning, all-in-one premier surgical platform	Our most popular and best value all-in-one platform	Excellence in rotary endodontics
	Recommended Use			
Oral Surgery		✓	✓	
Dental Implantology	✓	✓	✓	
Rotary Endodontics		✓	✓	✓
Restorative Dentistry		✓	✓	
	Motor and Handpiece Connection Type			
	Electric	Electric	Electric	Electric
	Special Features			
Multi-Function Foot Control	Yes or On/Off Option	Yes or On/Off Option	Yes or On/Off Option	
Handpiece Calibration	✓	Yes, Dynamometer	✓	✓
Handpiece Illumination		Yes, AEU-7000L-70V	✓	
Dentsply Sirona rotary endodontic file library		✓	✓	✓
User-customizable memory presets	✓	✓	✓	✓
Auto-stop-reverse for endodontics		✓	✓	✓
Upgradeable Software		✓	✓	✓
Autoclavable Tubing	✓	✓	✓	
	General Pricing*			
	\$3,300-\$3,500	\$5,230-\$5,750	\$3,730-\$3,950	\$2,300

CALL TODAY! 425.408.9347 or shop online at aseptico.com

1070 Series
Streamline & Simplify
Implant Procedures

AEU-1070-70V

The Aseptico 1070 Series Implant Motor is designed to simplify dental implant procedures by providing all the features most desired by today's implant dentists with a straightforward design, smaller footprint and a more affordable price.



Streamlined workflow

The 1070 Series Implant motor has the ability to place implants up to 70 Ncm with a 20:1 reduction handpiece and features a 1:1 setting for surgical and low-speed procedures.

User-friendly interface

The touchscreen display and intuitive user interface make it simple to adjust motor settings based on your preferences.

Hands-free control

The multi-function foot control enables easy adjustment of speed, direction, irrigation flow, torque, and presets, enhancing efficiency and infection control.

Precise and powerful

The brushless 40,000 RPM motor delivers the power needed for accurate and efficient implant procedures.

Versatility

With the 1:1 setting (and available 1:1 straight handpiece attachment), the motor is also suitable for surgical and low-speed procedures as well.

Customization options

The eight programmable presets allow you to customize the motor settings based on your preferred techniques.

Enhanced accuracy

The integrated handpiece calibration feature of the 1070 Series Implant Motor ensures speed and torque accuracy at the bur, providing consistent and reliable performance.

Compact and durable

The motor's compact footprint and water-resistant chassis design make it suitable for various clinical environments.

AEU-1070-70V

- Recommended for dental implants
- Up to 70 Ncm to streamline implant workflow
- Compatible with E-type 20:1 and 1:1 handpieces with external irrigation connections
- Advanced features and efficient adjustable torque desired by implant dentists
- Autoclavable 40k rpm brushless micromotor
- Works for all implant types
- Adjustable irrigation flow rate of up to 140 ml/min.
- Reusable, autoclavable irrigation tubing
- Easy cleaning and maintenance
- Console dimensions: 19.6cm W x 18.0cm D x 16.0cm H (7.72" W x 7.09" D x 6.30" H)
- Universal voltage 100-240V

1070 Motor with On-Off Foot Switch	AEU-1070	\$3,330
1070 Motor with Multi-Function Foot Control	AEU-1070-70V	\$3,550

Options		
Micromotor w/ 10ft Cable	AE-230A-40EXT	\$1,200
Multi-Function Foot Control w/ 10ft Cable	AE-70V2-10	\$660
On/Off Foot Control w/ 10ft Cable	AE-7PM-10	\$210

Replacement Parts		
Micromotor w/ 6ft Cable	AE-230A-40	\$1,100
Multi-Function Foot Control	AE-70V2	\$610
On/Off Foot Control	AE-7PM	\$200
Irrigation Bag Hanger Rod	461541	\$48
Power Cord	840107	\$23
Handpiece Cradle	461561-01	\$11
Irrigation Tubing Set	AE-24	\$80
Irrigation Pump Tubing Section (qty: 10)	AE-24-PUMP	\$26
Tubing and Accessories see page 19-20		

Blog:
FAQ on the new 1070
Series Implant Motor



Watch a Video
Demonstration



Mont Blanc AHP-85MB-X shown (not included). See details and other compatible handpieces pages 12-18.

7000 Series Motor All-in-One Surgical Platform

Oral Surgery, Implant Placement,
Endodontics & Restorative Dentistry

AEU-7000



AEU-7000L with Mont Blanc AHP-85MBFO-CX implant handpiece. Handpiece not included.
See details and other compatible handpieces pages 12-18.

The Aseptico 7000 Series Dental Motor

is one of the most powerful and easy-to-use implant/surgical motors available today. This all-in-one platform provides specialist-level performance for restorative, oral surgery, implant, and rotary endo applications.

All-in-One Platform

The 7000 Series has the ability to place implants up to 80 Ncm with a Mont Blanc® 20:1 handpiece.

For oral surgery, predictable and adjustable torque adds to practice efficiency.

For endodontics, the 7000 Series features auto-stop-reverse functionality to prevent file separation and comes preloaded with the complete Dentsply Sirona rotary file library plus the ability to customize your own presets.

For restorative procedures, the 7000 Series streamlines workflow; simply change out the handpiece and you have a restorative powerhouse.

Dynamometer Handpiece Calibration

The 7000 Series dynamometer calibration system gives dentists the ability to quickly and easily diagnose the health of a handpiece and verify that the handpiece is producing the correct torque for the application.

Upgradeable Software

Upgradeable software means a fantastic return on investment while ensuring that your future needs are also met.

Autoclavable Tubing

Reusable, autoclavable irrigation tubing means less waste, no more buying or keeping track of costly disposables.

Multi-Function Foot Control

The Multi-Function Foot Control replaces the need to adjust the unit manually during the procedure providing the operator variable speed control, forward-reverse, irrigation flow rate, torque adjustment, and memory functions.

Available handpiece ratios

1:1, 1:2, 1:3, 1:5, 8:1, 20:1

AEU-7000

- Recommended for oral surgery, dental implants, rotary endodontics & restorative dentistry
- All-in-one platform streamlines workflow
- Compatible with E-type 8:1 endo, 20:1 implant, 1:1, 1:2 surgical, 1:5, 1:3 high speed handpieces with external irrigation connections
- Additional handpiece ratios available via dynamometer calibration 4:1, 16:1, 32:1
- Advanced features and efficient adjustable torque desired by specialists
- Dynamometer provides exceptional handpiece testing and calibration
- Autoclavable 40k rpm brushless micromotor
- Works for all implant types
- Endodontic presets and complete Dentsply Sirona rotary endodontic file library onboard
- Upgradeable software
- Adjustable irrigation flow rate of up to 140 ml/min.
- Reusable, autoclavable irrigation tubing
- Easy cleaning and maintenance
- Optional: fiber-optic illumination, multifunction foot control, carrying case
- Console dimensions 9.98" W x 9.42" L x 5.1" H (25 x 24 x 13 cm)
- Universal voltage 100-240V

7000 Motor with On/Off Foot Control	AEU-7000E	\$5,230
7000 Motor with Multi-Function Foot Control	AEU-7000E-70V	\$5,450

7000 Motor with LED Illumination and Multi-Function Foot Control	AEU-7000L-70V	\$5,750
--	---------------	---------

Options

Carrying Case	410205-HNDL	\$950
Micromotor with 10ft cable	AE-230M-40EXT	\$1,250
Fiber Optic Micromotor with 10ft cable*	AE-230L-40EXT	\$1,420
On/Off Foot Control	AE-7PM	\$200
On/Off Foot Control with 10ft cable	AE-7PM-10	\$220
Multi-Function Foot Control with 10ft cable	AE-70V2-10	\$660

Replacement Parts

Fiber Optic Micromotor*	AE-230L-40	\$1,300
Micromotor with 6ft cable	AE-230M-40	\$1,150
Multi-Function Foot Control	AE-70V2	\$610
Calibration Bur	461558	\$16
Irrigation Bag Hanger Rod	461541	\$48
Power Cord	840079	\$23
Handpiece Cradle	461561-01	\$11
Irrigation Tubing Set	AE-23	\$80
Pump Section Tubing	AE-23-PUMP	\$26

More Tubing and Accessories *see page 19-20*

*Requires AEU-7000L Console



"This is the best implant motor on the market. I've been using it for 8 years now and I've never had an issue. Very happy with this reliable equipment."

John Downey, DDS



Watch a Video
Demonstration



Blog: Aseptico® 7000 Series
Dental Motor is Named Top
25 Implant Product



PACKAGE DISCOUNTS AVAILABLE.
CALL TODAY! 425.408.9347

6000 Series Motor All-In-One Platform Our Best Value!

Osteotomy & Implant Placement,
Rotary Endo & Restorative

AEU-6000

- Recommended for implantology and rotary endodontics
- All-in-one platform streamlines workflow
- Compatible with E-type 8:1 endo, 20:1 implant, 1:1, 1:2 surgical, 1:5, 1:3 high speed handpieces with external irrigation connections
- Integrated Handpiece Calibration
- Autoclavable 40k rpm brushless micromotor
- Works for all implant types
- Endodontic presets
- Upgradeable software
- Adjustable irrigation flow rate of up to 140 ml/min.
- Reusable, autoclavable irrigation tubing.
- Easy cleaning and maintenance
- Optional Multi-Function Foot Control
- Small footprint:
8.6 W x 9.2" L x 4.8" H (22 x 23 x 12 cm)
- Universal voltage: 100-240V



AEU-6000 shown with On/Off Foot Pedal and Mont Blanc Handpiece (Not Included)

"I have really enjoyed my Aseptico AEU-6000. It meets all my surgery needs from placing implants to removing full bony impacted wisdom teeth. I have found that this unit has more features and better torque than my old W&H."

Heath Snell, DDS
General Dentist

The AEU-6000 is our best-selling all-in-one motor because of its reliability and affordability. This versatile machine has the ability to place implants up to 50 Ncm with a Mont Blanc 20:1 handpiece. It performs rotary endo with auto-stop-reverse, and it can also be used as a backup for your restorative and oral surgery systems. No other motor on the market can give you these benefits at such an affordable price!

With upgradeable software, the AEU-6000 unit will be a workhorse for many years to come. Autoclavable irrigation tubing means you don't have to buy a new set for each procedure. Five preset buttons can be programmed separately for both implants and endodontics, reducing your procedure time dramatically.

Implant Motor w/ Multi-Function Foot Control	AEU-6000-70V	\$3,950
Implant Motor w/ On-Off Foot Switch	AEU-6000	\$3,730
Options		
Carrying Case	410205-HNDL	\$950
Micromotor with 10ft cable	AE-230M-40EXT	\$1,250
On/Off Foot Control	AE-7PM	\$200
On/Off Foot Control with 10ft cable	AE-7PM-10	\$220
Mult-Function Foot Control with 10ft cable	AE-70V2-10	\$660
Replacement Parts		
Micromotor with 6ft cable	AE-230M-40	\$1,150
Multi-Function Foot Control	AE-70V2	\$610
Irrigation Bag Hanger Rod	461541	\$48
Power Cord	840079	\$23
Handpiece Cradle	461561-01	\$11
Irrigation Tubing Set	AE-23	\$80
Pump Section Tubing	AE-23-PUMP	\$26

Tubing and Accessories see page 19-20



Electrical Operative Motor Top-Rated Motor Converts Air to Electric

AEU-5000



Recommended handpiece 1:5 Restorative with fiber-optic illumination (AHP-72MK-FO)

The versatile and cost-effective operative motor is designed to convert a standard air dental delivery system into a powerful electric handpiece system for numerous dental applications.

The AEU-5000 features endo functionality, fiber optic intensity adjustment, a compact and easy-to-clean console, a large backlit display, intuitive user interface, and a 40,000 rpm brushless autoclavable micromotor. Innovative LED technology produces some of the brightest illumination you'll find.

- Converts air to electric
- Powerful, high torque electric motor
- Attach to your air-driven high speed connection and use your rheostat
- Electrical option for efficient and quiet restorative procedures
- Up to 200,000 rpm with high speed handpiece
- Intuitive design and operation
- Includes endo mode with auto-stop-reverse torque setting
- Universal voltage 100-240V

Electrical Operative Motor	AEU-5000	\$1,890
Adjustable Mounting Bracket	330569	\$110
Replacement Micromotor	AE-240-40	\$1,550
1:5 Restorative Handpiece with Fiber-Optics	AHP-72MK-FO	\$990

Recommended Handpieces see pages 12-18

Available handpiece ratios 1:5, 1:2, 1:1, 5:1, 8:1



Optional AMC-25 Integration
Read more about the benefits of our most advanced dental cart on page 38.

AMC-25 Product Shown Not Included, Sold Separately

Endo DTC® Rotary Motor America's Most Popular

AEU-25

- Most sold, most trusted, endo motor in America
- Dentsply Sirona file library included
- Auto-stop-reverse functionality



Foot control included



MicroNiTi Handpiece

We recommend the AHP-88MNP MicroNiTi endodontic handpiece. See page 15.

The Endo DTC Rotary Endo Motor is America's most popular endo motor because it is a proven reliable performer. Many dentists did their endodontic training on this system and they still prefer it today. There is so much packed into this small footprint at such an affordable price: advanced digital torque control technology with auto-stop-reverse helps prevent file separation, five programmable memory buttons, onboard file library, integrated calibration system, autoclavable 30k rpm brushless micromotor, and upgradeable software. Pair with a AHP-88MNP handpiece to make the perfect rotary endo system.

Available handpiece ratios
1:5, 1:1, 4:1, 8:1, 16:1, 64:1

Endo DTC®	AEU-25	\$2,200
Replacement Parts		
Micromotor	AE-4B-30	\$790
On/Off Foot Control	AE-7P	\$190
Handpiece Holder	AE-5FG	\$25
Motor Cap (For AE-4B-10/30)	461084-08	\$10
Recommended handpieces see page 15		

Chad Burgess, DDS
Endodontist

"As a busy endodontist, I want equipment that is reliable, but more importantly, from a company that has great customer service. Aseptico delivers on both fronts."

Electronic Endodontic System Our Most Advanced Endo Motor

AEU-27A

- Lightweight motor and cable
- Upgradeable software
- Weighted omnidirectional foot pedal
- Includes all features of America's most popular endo motor, the Endo DTC®



The Electronic Endodontic System features digital torque control with auto-stop-reverse to prevent file separation and automatic calibration to test both the motor and handpiece for accurate performance and consistent operation. Easy to operate, with an auto-clavable brushless motor, this motor also features onboard file library, five programmable memory locations, and preset options that allow single file, file series, or custom file series. Upgradeable software assures return on investment. Includes weighted omnidirectional foot pedal.

Available handpiece ratios
1:5, 1:1, 8:1, 16:1

Electronic Endodontic System	AEU-27A	\$2,300
Replacement Parts		
Micromotor	AE-250-30	\$670
On/Off Foot Switch	AE-7PM	\$200
Motor Cap (For AE-4B-10/30)	461084-08	\$10
MicroNiTi endodontic handpiece recommended	see page 15	

Richard Mounce, DDS

"The AEU-27A is the most versatile, safe and economical motor available anywhere in the world, at any price. Upgradeable, simple to use, and ergonomic, it sets the standard for electric endodontic motors for many years to come."



Don't forget the HandiDam®
Created by Aseptico, it's the original pre-framed rubber dam trusted worldwide. Easier and more efficient — it's the only rubber dam that allows access to suction or take x-rays while isolation and clamps remain in place. See page 27.



Mont Blanc® World Class Implant Handpieces

The Mont Blanc Implant Handpiece is crafted from stainless steel and titanium, making it lightweight and durable. Precision ball bearing construction reduces friction to prevent overheating during surgical procedures. It features smooth surfaces that improve infection control and ensure easy cleaning, with the ability to remove the head for cleaning without tools.

Electric (E-type) | 80 Ncm max torque | 20:1 reducer | 15-2,000 rpm (40k rpm micromotor)



Economy Best Value Implant Handpiece

With an upper torque limit of 70 Ncm, the Economy Dental Implant Handpiece is ideally suited as a primary or backup handpiece for the AEU-1070 and AEU-6000 platforms. This implant handpiece offers exceptional value and is excellent for drilling osteotomies, tapping and threading.

Electric (E-type) | 70 Ncm max torque | 20:1 reducer | 15-2,000 rpm



AseptiSpray™ Ensuring Handpiece Longevity



Enhance Your Handpiece Performance and Longevity with AseptiSpray™ Advanced Cleaner and Lubricant:

- Provides complete internal care for all types of air-driven (pneumatic) and E-type (electric) handpieces.
- Made by Aseptico to ensure the highest standards for modern dental equipment.
- Features a CFC-free and ozone-safe formula for responsible practice.
- Proper lubrication enhances performance and extends the life of your handpieces.

Learn More On Page 19.

Watch a Video
Demonstration



ASEPTICO SELECT HANDPIECES ARE EXCEPTIONAL PRODUCTS,
RIGOROUSLY QUALIFIED AND PERFORMANCE TESTED.

See all handpiece accessories, including AseptiSpray™ handpiece cleaner and lubricant, irrigation clips, cannulas, and more on pages 19-20.

1:1 Straight Surgical

AHP-64MKI

- **Internally Irrigated**
- Electric (E-type)
- 300 to 40,000 rpm (w/ 40k rpm micromotor)
- 1:1 low speed



MK-dent 1:1 Straight Electric Surgical Handpiece is ideal for osteoplasty, working dentin, adjusting dentures, sectioning enamel and performing sinus lifts. The spray nozzle ensures safe working conditions due to the optimum cooling of the preparation site. The nozzles are made to provide a clear view of the surgical area as well as provide a pleasant noise level.

Straight Electric Handpiece	AHP-64MKI	\$650
AseptiSpray™ Cleaner and Lubricant	AHP-09	\$44

20° Surgical

AHP-77W

- Electric (E-type)
- 20° bend for optimum comfort
- 1:2 increaser
- 600 to 80,000 rpm (w/ 40k rpm micromotor)
- Requires 70 mm bur length



AHP-77W 20° Electric High Speed Handpiece is the oral surgeon's number one choice for third molar extractions. The angled design prevents fatigue by allowing it to rest on the back of the hand. A scratch-resistant surface and reinforced nose bearings guarantee durability, and the sliding lip seal protects the interior from penetration of dirt particles.

20° Electric Handpiece	AHP-77W	\$2,047
W&H Cleaner and Lubricant	AHP-10 & AHP-10N	\$76.75
Irrigation Clip	AHP-77W-CLP	\$69.75

45° Surgical

AHP-79W

- Electric (E-type)
- 1:2.7 increaser
- 810 to 108,000 rpm (w/ 40k rpm micromotor)



W&H 45° High Speed Surgical Handpiece is ideal for precise tooth separations and apical resections. The 45° head allows for optimum access to the posterior of the mouth. Available in non-fiber optic and self-generating fiber-optic options, this handpiece comes complete with a triple-spray head that guarantees sufficient cooling of the surgical bur and site.

45° Electric Handpiece	AHP-79W-FO	\$2,123
w/LED		
45° Electric Handpiece	AHP-79W	\$1,761
W&H Cleaner and Lubricant	AHP-10 & AHP-10N	\$76.75

1:2 Straight Surgical

AHP-77DS

- Electric (E-type)
- 1:2 increaser
- Up to 80,000 rpm



AHP-77DS Electric High Speed Handpiece is lightweight and compact. The traditional straight design allows for exceptional posterior access. With its 1:2 increaser and rotation speeds up to 80,000 rpm, this handpiece is a great option for bone contouring and sectioning teeth.

1:2	AHP-77DS	\$1,950
Replacement Irrigation Clip	AHP-77D-NOZ	\$32
AseptiSpray™ Cleaner and Lubricant	AHP-09	\$44

Microhead Endodontic

AHP-88MNP

- Electric (E-type)
- Push button
- Speed range (w/ 30k rpm micromotor)
- 8:1 — 40 to 3,600 rpm



Anthogyr MicroNiTi Microhead Endodontic Handpiece has a push button head and provides the same unrivaled accuracy as other MicroNiTi series handpieces, whether it is being used with a rotary or reciprocating motor. The micro-head size allows for exceptional access and line of sight during procedures. Made for nickel-titanium endodontic files, the MicroNiTi Microhead is compatible with all ISO rotary endo files.

8:1	AHP-88MNP	\$780
AseptiSpray™ Cleaner and Lubricant	AHP-09	\$44

1:3 Restorative

AHP-65TI-FO and AHP-65TI

- Electric (E-type)
- External irrigation
- 1:3 increaser
- Optional fiber optics
- Max 120,000 rpm (w/ 40k rpm micromotor)
- Perfect for high speed dentistry with your implant/surgical motor



With an external irrigation port and 1:3 increaser gearing, the AHP-65TI Restorative Electric High Speed Handpiece makes any surgical motor with irrigation into a restorative powerhouse. Ergonomic and light-weight design allows for optimum comfort during your procedure.

1:3 High Speed w/ Fiber Optics	AHP-65TI-FO	\$1,700
1:3 High Speed	AHP-65TI	\$1,150
AseptiSpray™ Cleaner and Lubricant	AHP-09	\$44

1:5 Restorative

AHP-72MK-FO

- * Electric (E-type)
- Internal irrigation
- 1:5 increaser
- Fiber optics
- Max 200,000 rpm (w/ 40k rpm micromotor)
- Perfect for high speed dentistry with your electric operative motor



Provides powerful, quiet, vibration-free operation with class-leading precision, the 1:5 Restorative Handpiece features Long Lasting Gears that are engineered to excel at every angle, and ceramic bearings that stand up to heavy use. Fiber optics optimize visibility with 100% glare-free illumination at a color temperature of 5500 K.

1:5 High Speed	AHP-72MK-FO	\$990
AseptiSpray™ Cleaner and Lubricant	AHP-09	\$44

4-Hole Restorative

AHP-96MK

- Air-driven
- Internal irrigation
- Non-optic
- Max 320k rpm, 24 watts
- Max volume 62 dB
- Recommended for high speed dentistry



Simple to use technology, developed for what you do every day. The 4-hole High Speed Handpiece features a clean, elegant shape and a virtually scratch proof smooth chrome surface. Manufactured using proprietary laser balancing technology, this handpiece delivers 24 watts of torque and features a premium cartridge for enduring precision and prolonged life.

MK Dent 4-hole High Speed Restorative Handpiece	AHP-96MK	\$500
AseptiSpray™ Cleaner and Lubricant	AHP-09	\$44

5-Hole Restorative

AHP-97MK

- Air-driven
- Internal irrigation
- Fiber optic glare-free, color accurate illumination
- Max 320k rpm, 24 watts
- Max volume 62 dB
- Perfect for high speed dentistry



It feels great in your hand and sports a clean modern shape with a chrome surface. Fiber optics deliver 100% glare-free, color accurate illumination. The 5-hole High Speed Restorative Handpiece is equipped with the AeroPerformance Blade premium cartridge. Ceramic bearings stand up to heavy use and autoclaving.

MK Dent 5-hole High Speed Restorative Handpiece	AHP-97MK	\$530
AseptiSpray™ Cleaner and Lubricant	AHP-09	\$44

6-Pin Restorative

AHP-82MK

- Air-driven
- Internal irrigation
- Fiber optic with LED color accurate illumination
- Max 350k rpm, 25 watts
- Max volume 55 dB



Excellent for high speed dentistry in-office or with a mobile pneumatic delivery system. More power, better sound, spectacular control. The MK Dent 6-pin High Speed Restorative Handpiece is built for the pro user who demands best in class. Reverse flow protection minimizes the entry of contaminants into the head. Fiber optics provide color-accurate illumination. Supremely finished with a scratch-resistant titanium surface. Dentists will enjoy working with this exceptional handpiece.

MK Dent 6-pin High Speed Restorative Handpiece	AHP-82MK	\$1,050
MK Dent 6-pin Coupler	AHP-82MK-QD	\$390
AseptiSpray™ Cleaner and Lubricant	AHP-09	\$44

6-Pin Fiber Optic Coupler

AHP-82MK-QD

- Air-driven
- Works with LED bulbs and glass rod fiber optics
- Color accurate



AHP-82MK Handpiece shown with AHP-82-MK-QD Coupler attached.

This 6-Pin Fiber Optic Coupler comes with a non-retraction valve and water volume regulation. Works with both LED and glass rod fiber optics to provide 100% glare free and color accurate light.

6-Pin Fiber Optic Coupler AHP-82MK-QD \$390

CALL TODAY!
425.408.9347
Or shop online at
aseptico.com

Low Speed Straight

- AHP-101
- Electric (E-type)
 - Attaches to E-type air motor
 - Also pairs with AHP-50MK
 - 1:1 Straight Handpiece
 - Post for Doriot-style prophy angle attachment
 - Recommended for use with Aseptico portable systems and lab motors



The Low Speed Straight Handpiece is made of lightweight aluminum and is designed to work with Aseptico® air-driven portable dental systems. The AHP-50 Low Speed Air Motor, which allows air turbines to accept E-type handpieces, is required to connect this handpiece to Aseptico portable dental systems.

Low Speed Straight Handpiece **AHP-101** \$380

Low Speed Contra Angle

- AHP-63MK
- Electric (E-type)
 - Attaches to E-type air motor
 - Also pairs with AHP-50MK
 - Swing Latch
 - 1:1 Contra Angle
 - Recommended for low speed work on portable equipment or in an office operatory



Low Speed Contra Angle Handpiece is a D-latch style low speed handpiece that can be attached to any E-type air motor. Made by MK Dent

MK Dent 1:1 Swing Latch Contra Angle	AHP-63MK	\$270
AseptiSpray™ Cleaner and Lubricant	AHP-09	\$44

4-Hole to E-Type Air Motor

- AHP-50MK
- Air Driven
 - Low Speed Air Motor (4-hole)
 - Adapts an air system to accept an E-type handpiece



The Low Speed Air Motor accepts E-type handpieces. Recommended for use with Aseptico® air-driven portable dental systems or other in-office applications where you need to adapt an air system to an E-type handpiece.

Low Speed Air Motor **AHP-50MK** \$380



AseptiSpray™ Handpiece Cleaner and Lubricant

AseptiSpray™ is a must for proper cleaning and maintenance of handpieces. It keeps handpieces working efficiently and ensures long life and durability.

AseptiSpray™ w/ Nozzle (8.8 oz can)	AHP-09	\$44
Nozzle only	AHP-09N	\$15

W&H f1 Service Oil Handpiece Cleaner and Lubricant

Synthetic oil of maximum purity, which has been specially made for maintenance of W&H handpieces. Required by W&H warranty policy

Handpiece Service Oil (Spray Cap NOT included)	AHP-10	\$65
Service Oil Spray Cap	AHP-10N	\$11.75

Blog: Advantages of Autoclavable Tubing for Oral Surgery Systems



Mont Blanc®/Economy Handpiece Accessories

Depth Stop Attachment	AHP-85MB-STP	\$130
Black Irrigation Clip	AHP-85MB-CLP	\$55
Cannula	AHP-85MB-CAN	\$30
Cannula Clip Set (Includes Small Clip, Large Clip, Cannula and Tubing)	AHP-07	\$24
Cannula and Tubing	AHP-04	\$19
Small Clip	AHP-05	\$8
Large Clip	AHP-06	\$8
Y Connector	AE-2	\$9



Autoclavable Black Irrigation Clip for Mont Blanc® Handpiece
AHP-85MB-CLP

Tubing & Accessories

Autoclavable surgical tubing saves money and prevents waste. Disposable tubing sets cost about \$12 per set compared to less than \$1 per use for autoclavable tubing that lasts for approximately 100 uses. Autoclavable irrigation clips are money-savers too.

Irrigation Tubing Set for IV Bags (For AEU-6000 and 7000 series)	AE-23	\$80
Irrigation Pump Tubing Section (qty:10) (For AEU-6000 and 7000 series)	AE-23-PUMP	\$26
Irrigation Clips (qty: 10)	AHP-03	\$22
Irrigation Pump Tubing Section (qty: 10)	AE-23-10	\$110
Irrigation Tubing Set for IV Bags (For AEU-1000 and 1070 series)	AE-24	\$80
Irrigation Pump Tubing Section (qty:10) (For AEU-1000 and 1070 series)	AE-24-PUMP	\$26
Motor Cap (For AE-4B-10/30)	461084-08	\$10
Motor Plug (For AE-4B-10/30)	460183-08	\$11

CALL TODAY! 425.408.9347 or shop online at aseptico.com

Foot Controls

Multi-Function Foot Control

AE-70V2

The Multi-Function Foot Control replaces the need to adjust the unit manually during procedures, providing the operator with variable speed control or on/off pedal, forward-reverse, torque adjustment, irrigation flow rate, and memory functions for implants and endodontics. The Multi-Function Foot Control is available with a standard 6 ft. cable (AE-70V2) or a 10 ft. cable (AE-70V2-10). Made in USA.

Multi-Function Foot Control with 6' cable	AE-70V2	\$610
Multi-Function Foot Control with 10' cable	AE-70V2-10	\$660

On/Off Foot Switch



AE-7PM

On/Off Foot Switch for Aseptico dental motors:
AEU-7000 Series, AEU-6000 Series, AEU-27A, and portable systems AEU-350, AEU-525/AEU-525S. Made in USA.

Foot Switch with 6' cable	AE-7PM	\$200
Foot Switch with 10' cable	AE-7PM-10	\$220

Replacement Micromotors

AEU-7000L micromotor and 6' fiber optic cable	AE-230L-40	\$1,300
10 ft extended cable option	AE-230L-40EXT	\$1,420
6000 Series and 7000 Series Implant Motors micromotor and 6' cable	AE-230M-40	\$1,150
10 ft extended cable option	AE-230M-40EXT	\$1,250
AEU-27A Advanced Endodontic Motor micromotor	AE-250-30	\$670
AEU-5000 Electric Operative Motor micromotor	AE-240-40	\$1,550
Micromotor and Cable for:	AE-4B-30	\$790
AEU-25 Endo DTC®		
AEU-707AV2 Implant/Surgical Motor		
AEU-17BV2 SurgiMotor II		
Motor Cap for AE-4B-10/30	461084-08	\$10
AEU-10SS Light Duty Lab Motor replacement micromotor	AM-103B	\$220
AEU-1000/1070 Micromotor	AE-230A-40	\$1,100
AEU-1000/1070 Micromotor with 10 ft Cable	AE-230A-40EXT	\$1,200



Aseptico
INNOVATIONS
AND UTILITY





Newest Generation Penguin II

Implant Stability Monitor

RFA-55250-MULTI

New and Better! Aseptico is proud to announce the US arrival of the next generation RFA implant stability monitor – the Penguin II.

- Evolved design brings increased comfort, surer grip
- Relocated, easier-to-read LED display screen
- Battery can now be replaced, extending product life
- New vertical wireless charging stand saves space

Spotter Smart Implant Detector Find Buried Implants Faster and Preserve Soft Tissue

SPT-5I11001-KIT

- Painlessly and accurately locate buried implants quickly and easily
- Non-invasive, no exploratory flap needed
- Minimize the recovery time for the patient
- Two-year warranty
- Works with any titanium implant


Spotter Smart Implant Detector locates dental implants covered by gum tissue quickly and painlessly — no exploratory flap needed! Just move the wand over the site and the Spotter will guide you to the exact center of the implant.

Instead of opening a flap to locate the implant, the Spotter finds the exact center quickly, allowing you to punch and reveal only the implant, thereby helping to preserve soft tissue and minimize recovery time for the patient.

The Spotter is especially helpful in the case of multiple implants when implant sites may be hard to visualize. The Spotter Smart Implant Detector works with any type of titanium implant.

It is easy to use, requiring only a few minutes to learn how to operate. The unit operates on one AA battery and its detachable sensing head is autoclavable.

Spotter SPT-5I11001-KIT \$700
Sensors 5/PKG SPT-5I11101 \$130



"With Spotter, we don't have to raise a flap, we can go in a flapless procedure, just punching the soft tissue to get to the center of the implant."

Ziv Mazor, DMD

spotter
 Smart Implant Detector

Implant Stability Monitors

- Provides accurate Implant Stability Quotient measurement for greater confidence determining when to load dental implants
- Cordless, lightweight, and simple to use
- Half the price of competing ISQ devices
- Two-year warranty, Made in Sweden

Penguin Implant Stability Monitors

Now every dentist can access affordable Resonance Frequency Analysis (RFA) technology with the Penguin RFA and Penguin II for assessing osseointegration of dental implants. These devices measure local bone density, translating it into an implant stability quotient (ISQ) value between 1 and 99. Higher ISQ values indicate greater stability. Especially useful for all-on-4 or all-on-X cases, the Penguin monitors provide precise, objective insights into osseointegration, ensuring predictable outcomes.

MulTipeg™ Sensing Devices

Both Penguin RFA and Penguin II utilize the same MulTipecs™

- For all major implant systems
- Made of durable, tissue friendly titanium
- Autoclavable for 20 uses
- Laser marked with type number
- ISQ Standard Calibrated™
- Cost effective

Penguin II	RFA-55250-MULTI	\$2,490
Penguin RFA	RFA-55002-MULTI	\$2,290
Penguin MulTipeg	RFA-550**	\$39
MulTipeg Driver	RFA-55003	\$25
Sterile Covers (20 Pack)	RFA-55105	\$44
Penguin II Replacement Battery	RFA-55291	\$200

**Part number varies by implant platform.

Watch a Video
 Demonstration



View MulTipeg List



CALL TODAY! 425.408.9347 or shop online at aseptico.com

Made in USA Trusted Worldwide.

In 2025 Aseptico Inc. reaches the milestone of 50 years producing innovative dental equipment for dental professionals around the world. As a company, we are privileged to serve the dental profession and value Aseptico's significant role in the ongoing advancements of rotary instrumentation — especially in the fields of endodontics and dental implantology. In our role as world leader in mobile & portable dental equipment, Aseptico is passionate about innovating ways to help dental care reach more people in more places than ever before. We love our job solving critical equipment challenges to advance dental care!

Almost Everything Has Changed

The world has changed dramatically since Aseptico was founded in 1975. The first Apple computer was introduced in 1976. The dental world saw breakthroughs around the same time. In the late 70's, Dr. PI Brånemark's work accelerated the advancement of dental implants, setting the stage for decades of innovation that is still charging forward.

In the 1990s, the introduction of engine-driven, nickel-titanium (NiTi) rotary files for endodontic procedures transformed root canal therapy.



Aseptico was, and still is, at the forefront of designing and engineering the rotary instrumentation dentists need to perform root canals and get patients out of pain.

It's a Mobile World!

When Aseptico introduced its first portable dental system, a state-of-the-art portable system was made up of 66 pounds of equipment in two boxes. Today, our GO Ultra Portable dental system is so light and compact that it qualifies as carry-on luggage; it features an electric motor with LED illumination, premium ultrasonic scaler, and can operate on electricity anywhere in the world — or if there is no electrical grid, it can operate on external 24V battery power.

50
YEARS
OF EXCELLENCE
ASEPTICO

Commitment to Excellence, Patient Care

Aseptico is committed to excellence in all that we do and every product we make because dental clinicians and their patients are relying on us every day. Accordingly, Aseptico® products are fully regulatory compliant to US and international manufacturing and safety standards, including FDA Good Manufacturing Practices, ISO 60601-1, ISO 13485:2016, UTL/ETL certification for electrical, mechanical, and fire safety, and CE product certifications.

Trusted Worldwide, Forward Leading

Exciting new advancements are taking place in dentistry every day that keep us innovating — and we wouldn't have it any other way. The Aseptico team is poised to engineer and deliver the state-of-the-art dental equipment needed to advance clinical outcomes around the world.

Ready for the Future!

As Aseptico celebrates this 50-year milestone, its commitment to the dental profession is stronger than ever. From cutting-edge dental motor systems for oral surgery, implantology, endodontics and restorative dentistry, to the world's most advanced portable dental systems, and all the other innovative products in its Global Dental Catalog, Aseptico's reputation for innovation, quality and reliability is unmatched. Trusted worldwide, proudly made in USA, and driven by a passion to serve dental professionals, Aseptico is ready for the future of dentistry.

2024 PRODUCT RECOGNITION

Jan/Feb	AMC-15 Mobile Dental Cart	TOP INNOVATIVE PRODUCT
March	7000 Series Motor	TOP 25 IMPLANT PRODUCT
April	7000 Series Motor	TOP 25 ENDODONTIC PRODUCT
July/Aug	1070 Series Implant Motor	TOP 100 PRODUCT
Oct	7000 Series Motor	EDITOR'S CHOICE AWARDS 2024
Nov/Dec	GO Ultra Portable Dental System	TOP 50 TECHNOLOGY PRODUCT





Powered Trolleys

ATC-12V2 and ATC-15

- Keeps equipment organized and solves operatory space concerns
- Storage bin for foot pedals or other items
- Dual outlet with coiled AC power cord
- Adjustable height to easily store under a counter: 26.4" to 37.4" (67-95 cm)
- Four shelves w/ 3 removable stainless steel trays

Aseptico Powered Trolleys give dental teams the ability to easily roll dental motors, scalers, lasers, surgical systems or other electrical equipment from operatory to operatory. Available in large and standard sizes, each Powered Trolley is adjustable in height and features four shelves of storage with removable stainless-steel trays that fit securely into each shelf. High-quality, sturdy trolleys move and turn with ease and fit into almost any small space.

Standard Powered Trolley **ATC-12V2** \$795

Cart size: 26.4-37.4" x 15.75" x 16.5" (67-95 x 40 x 42 cm)
Tray size: 11.3" x 9" (28.8 x 22.8 cm)

Large Powered Trolley **ATC-15** \$950

Cart size: 26.5-38" x 19.5" x 19.75" (68-97 x 50 x 50 cm)
Large Tray: 15.4" x 11" (39.1 x 28 cm)

Powered Aseptico Trolleys

Featuring two heavy duty electrical outlets that are certified to UTL/ETL safety standards



Deluxe Operatory Cart

ATC-16

Revolutionize efficiency and organization in your dental practice with the Deluxe Operatory Cart. This compact, mobile cart features durable construction and smooth-rolling casters for effortless movement. With ample storage—multiple drawers and adjustable shelves—its ergonomic design ensures quick access to instruments, optimizing workflow and maximizing operatory space.

Deluxe Operatory Cart **ATC-16** \$4,295

Cart size: 16.5" W x 18.75" L x 34.75" H

Accessories


Adjustable, Articulating Arm	AA-16-ARM	\$900
Motor and IV Bag Holder	AA-16-IMPACC	\$750
Piezo Handpiece Holder	AA-16-PZACC	\$370

Products Shown Not Included, Sold Separately

HandiDam® pre-framed dental dam



- Pre-framed, flexible design facilitates access to the oral cavity for suction or X-rays
- Superior elasticity and tear-resistance
- Saves time—eliminates the need to remove and replace traditional dam during the procedure
- Ergonomic and comfortable for greater patient acceptance
- Meets FDA Powder-Free standard
- Extremely low protein content reduces patient irritation (<50 micrograms)
- Available in Latex or Latex-Free



HandiDam® Latex (Blue)		
Qty 20,	I-70BM	\$43
Qty 100,	I-70BM-100	\$165
HandiDam® Latex-Free (Mauve)		
Qty 20,	I-70LF	\$68
Qty 100,	I-70LF-100	\$280


TECHNIQUE CHAIN & GRID



Technique Chain & Grid easily attaches to a pre-emerged tooth, so that an orthodontic device may be used to draw the tooth into proper position. Aseptico is the originator of this remarkable traction system that professionals use worldwide.

- 10 sets per box
- 14k gold chain
- 3 grid shapes available

Cuspid (qty: 10)	AMA-09A	\$365
Bicuspid (qty: 10)	AMA-09A2	\$365
Incisor (qty: 10)	AMA-09A3	\$365



Dental Dam Clamps

Our five most popular clamps.
See full selection at aseptico.com
Dental Dam Clamp \$17 each



K-211
Universal labial cavities on anteriors



K-5
For lower molars



K-9
Most popular anterior clamps



K-90N
Most popular anterior clamps



K-8A
For lower molars, with strong festooned beaks





Automated Air Delivery System

ADU-07 and ADU-07A

- Choose 2 or 3 handpiece system
 - Includes 3-way air/water syringe
 - Wall or cabinet mount
- Optional fiber-optics
 - Includes foot control

The Automated Air Dental Delivery System is a pneumatic handpiece control system that gives you everything you need at a price you can afford. Highly customizable, select from a 2-handpiece system with manual control, or 3-handpiece system with automatic selection and individual air and water adjustments. Both models include air pressure gauge, 3-way air/water syringe, foot control, and water flow adjustments. Flat top design is ideal for mounting into cabinets. Requires connection to air and water sources.

Options

2-Handpiece/Manual Control	ADU-07	\$1,100
3-Handpiece/Automatic	ADU-07A	\$1,200
5-hole Fiber Optic Upgrade	AA-19LED-04	\$720
6-pin Fiber Optic Upgrade	AA-19A-06	\$720

Replacement Parts

Disc Foot Control w/ Water Toggle	AA-43W	\$110
Air/Water Syringe	TA-90D	\$85
Autoclavable Tips for Air/Water Syringe	TA-1	\$9



**Folding Arm
Wall Mount**

Folding wall mount for ADU-07,
ADU-07A and ADU-04, Dental
Delivery Systems

AA-05A \$450



Manual Air Delivery System

ADU-04

- 2-handpiece control system
 - 3-way air/water syringe
 - Optional independent water
- Ideal for arm or wall mounting
 - Requires connection to air

The Manual Air Delivery System is a pneumatic handpiece control system for two handpieces. It includes a 3-way air/water syringe and disk foot control. Requires connection to air and water sources or use optional independent water.

Manual Air Control System	ADU-04	\$650
Clean Water System with 1 Liter Bottle (pictured)	NWS-71	\$130

Replacement Parts

Disc Foot Control	AA-42	\$65
Air/Water Syringe	TA-90D	\$85
1L Water Bottle (for NWS-71)	730471	\$13
Autoclavable Tips for Air/Water Syringe	TA-1	\$9



MOBILE & PORTABLE DENTAL EQUIPMENT



FRIENDLY AND PROFESSIONAL SERVICE. CALL TODAY! 425.408.9347



Dental Care Anywhere

Aseptico® Mobile & Portable Solutions

Mobile Dental Carts

Portable Delivery Systems

Patient & Clinician Seating

Portable X-ray

Portable Exam Lights



NEW
AEU-350S

Featuring
**ULTRASONIC
SCALER**



Portable Seating Up to Any Task!

Aseptico®
Trusted Worldwide

MOBILE & PORTABLE SOLUTIONS Dental clinicians, public health outreach programs, humanitarian dental groups, and the US military rely on Aseptico® mobile and portable dental equipment to deliver dental care around the world. **Made in USA.**

FRIENDLY & PROFESSIONAL CUSTOMER SERVICE CALL TODAY! 425.408.9347
or shop online at aseptico.com



Aseptico®
Trusted Worldwide

Find The Right Portable Equipment For You
Driven to advancing dental care and patient access to care. **Call us today! 425-408-9347.**



	GO Ultra Portable	Express Portable	Transport III Portable	Lightweight Portable	AMC-15 Mobile Dental Cart	AMC-25 Mobile Dental Cart
	Our most advanced technology	World's #1 selling portable system	Military-grade performance	Only 17 pounds! Requires Portable Compressor	Affordable workhorse for general dentistry	Highly customizable, Hospital-grade dental cart
	Recommended Use					
General Dentistry	✓	✓	✓	✓	✓	✓
Dental Hygiene	✓	✓	✓	✓	✓	✓
Dental Sealant Placement	✓	✓	✓	✓	✓	✓
Rotary Endodontics			✓			Optional
	Motor and Handpiece Connection Type					
	Electric	Air	Electric	Air	Air	Air with Electric Option
	Special Features					
Handpiece Illumination	✓	Optional	✓	Optional	✓	✓
Integrated Ultrasonic Scaler Option	AEU-350S		AEU-525S			✓
Integrated Curing Light Option						✓
	Weight					
	32 lbs.	43.5 lbs.	56 lbs.	17 lbs.	139 lbs.	167 lbs.
	Qualifies as a carry-on luggage	Can be checked as luggage	Extra luggage fee required	Check with your airline	Usually moved by truck or van then rolled room-to-room in institutional settings	
	General Pricing*					
	\$7,850-\$9,500	\$5,850	\$11,700-\$13,500	\$1,990+	\$12,500	\$15,600



Perfect with the NEW GO Cart
Learn more about the GO cart on page 34.

NEW
ULTRASONIC
SCALER OPTION



Designed for Ease-of-Use and Reliability

The GO leverages Aseptico’s expertise in crafting leading portable dental units for today’s challenges and opportunities. Designed for operator ease-of-use and outstanding reliability, the GO is capable of restorative dentistry and hygiene functions.

At only 32 pounds, the GO Ultra Portable dental system is so compact it qualifies as carry-on luggage for air travel.

The GO delivers powerful airflow and suction yet operates more quietly and cleanly than standard air compressor units (it even has an air filter to reduce airborne particulates).

GO Ultra Portable Dental Delivery System	AEU-350	\$7,850
NEW GO With Scaler	AEU-350S	\$9,500

Options

24 VDC Battery	AA-350-BATCAB2	\$110
Power Cord		

Replacement Parts For AEU-350S

Thin Subgingival	ASC-10-PE37	\$130
Slim Universal	ASC-10-PE38	\$130
Power Universal	ASC-10-PE39	\$130
Scaler Holder	AA-350-HLD1	\$34
Scaler Wrench	ASC-10-UNIV	\$55
Scaler Ergo Sleeve	ASC-10L-ERGOSLV	\$55

Replacement Parts

Micromotor with 5' cable	AE-242SC-40	\$1,350
Water Supply Bottle	AA-350-BOT	\$46
Waste Collection Tank	330811	\$95
Vacuum Hose Plug	AA-350-VPLUG	\$33
Vacuum Line Filter	AA-350-VLF	\$8
Air Filter	AA-350-FIL	\$37
HVE Valve and Hose Assembly	AA-350-HVE	\$160
Saliva Ejector Valve and Hose Assembly	AA-350-SAL	\$60
Foot Switch	AE-7PM	\$200
Power Cord	840041	\$20



Easy and Light
Only 32 pounds!

Quieter Operation

A quieter vacuum system that turns on only when HVE and saliva ejector are removed from their holders creates a quieter experience for both operator and patient.

500 mL Water Capacity

Plenty of water capacity for restorative cases

1 Liter Waste Capacity

Easy, clean removal of waste

Impact Resistant

Durable, wheeled case with comfortable handles

Weighted & Omnidirectional Foot Pedal

Stays in place and can be activated from any angle

Powerful Electric Motor

5k to 200k rpm range is ideal for prophyl and restorative.

Faster Drying for Sealants

- Syringe Performance:
- 80 psi initial
 - Up to 50 psi continuous

More Powerful Suction!

New vacuum blower technology

NEW Ultrasonic Scaler Option Now Available

Fully integrated piezoelectric ultrasonic scaler with LED illumination

Super Compact!

22" x 13" x 9.2" fits into the airline overhead as carry on luggage



Plug in Anywhere in the World

100-240 VAC & 24VDC input options

Can Operate on External Battery Power

Powering the GO from a 24 VDC battery is an option for those who work off-the-grid. The cable is purchased from Aseptico, and the customer supplies the 24V battery.



Watch Cascade Hygiene's Testimonial on the GO

"It's really user-friendly and quiet. Everything takes half the amount of time to do which allows me to spend more time with my patients. I've never worked with a company that cared so much about keeping me in business so that I can help take care of other people."

Erin Poicus, RDH
Cascade Hygiene



OUR MOST ADVANCED TECHNOLOGY!
CALL TODAY! 425.408.9347

Portable Supply Cart
GO Cart

ATC-350ZT

Move equipment effortlessly with the lightweight Aseptico GO Cart. Designed for mobile dental professionals, it offers smooth all-terrain mobility, organized storage with six pouches, and built-in water bottle holders. Weighing only 16.5 lbs, it also features a durable frame that doubles as a portable seat.

GO Cart	ATC-350ZT	\$705
Cart size: 24.5"H x 22"W x 14.5"D (including wheels) (51"H with fully extended handle)		

"The GO Cart has completely transformed my practice. It's incredibly durable, keeps me organized, and makes moving between patients a breeze,"

Sonya Dunbar, MPH, RDH
Geriatric Toothfairy | Mobile Dental Xpress



Durable Telescoping Handle
Extends Up to 51".

Perfect With The GO
Pair the GO Cart with the GO (AEU-350/s) and bring everything you need in one trip.

Extra Storage
Use the two additional side pockets for extra items.

Patented Steel Frame
Doubles as a portable seat, supporting up to 300 lbs for on-the-spot convenience.

All-Terrain Tires
Ensure smooth transportation over any terrain, making it ideal for remote and diverse settings.



Storage For Everything
Equipped with six storage pouches and customizable organizer boxes, perfect for storing handpieces, scalers, and other dental tools.



Two water bottle holders at the back can also carry sanitation supplies

Express Portable
Dental System

ADU-17A and ADU-17A2

ADU-17X and ADU-17X2

- Fast and efficient for dental sealants
- Available with one or two air handpiece connections
- High volume evacuator and saliva ejector
- Built-in maintenance free compressor
- Water port for Cavitron/scaler
- High capacity 1 gallon waste container
- 110V or 230V switchable voltage
- 44-48 lbs (20-21.8 kg)
- 10.5" x 14.5" x 24" (27 x 36 x 61 cm)

Aseptico's most popular portable dental systems, the Express Air and Express II portable dental systems are fully self-contained, lightweight, reliable performers that provide everything a dentist or hygienist needs to perform effectively in the field. These powerful units feature up to two air-driven handpiece connections, high volume and saliva ejector suction, internal water reservoir, 3-way syringe, and integrated air compressor. They are transported using a pullout handle and built-in wheel assembly. The system's compact size allows for convenient storage, while the impact-resistant case makes the Express ideal for any field environment. Designed for global use, all Express systems are compatible with 110V or 220V power sources and come in under the airline baggage weight limit.

Express Air – features variable rheostat foot pedal. The handpiece only operates when the variable rheostat foot pedal is depressed; handpiece operates independently from suction and 3-way syringe. **(Dentist preferred)**

Express Air one handpiece connection	ADU-17A	\$5,850
Express Air two handpiece connections	ADU-17A2	\$5,950

Express II – features on/off electric foot switch. All components are activated by the on/off foot switch. This configuration is sometimes desirable for special needs patients or other specialized applications. Please call for more information. **(Hygienist preferred)**

Express II one handpiece connection	ADU-17X	\$5,850
Express II two handpiece connections	ADU-17X2	\$5,950

Fiber Optic Upgrade Kits		
5-hole	AA-19LED-04	\$720
6-pin	AA-19A-06	\$720
6-pin LED handpieces	AA-19GEN-04T6*	\$260
*Will not function with incandescent bulbs. Consult your Aseptico sales rep for details.		

Handpieces sold separately
Recommended handpieces see pages 16-18.



"We love the ease of this piece of Aseptico equipment! It's essential for our business to run efficiently being able to utilize everything we would have in a private office setting while providing care in long-term care facilities. No sacrifice in power or suction either!"

Liz LaJoie, RDH
Director of Clinical Field Teams
Aria Care Partners

Transport III Portable Dental System

AEU-525S and AEU-525

- High torque electric motor with restorative and endodontic programs
- Fiber optic electric handpiece connection
- High volume evacuator and saliva ejector
- Quiet, intuitive, powerful
- Built-in maintenance free compressor
- Water port for Cavitron/scaler
- 115V or 230V Dual voltage
- Optional fiber optic piezo scaler
- Optional amalgam separator kit
- 57-63 lbs (25.8-28.6 kg)
- 17.9" x 10.4" x 22.1" (45.5 x 26.4 x 56.1 cm)


The Transport III is Aseptico's most capable portable dental system. Its features include noise reduction technology, an intuitive user interface, and 40k rpm brushless autoclavable micromotor. It is easy to set up, operate, and tear down. All components fit neatly into the impact resistant carrying case with inset handle and wheels.

Delivery module and user interface panel mount externally, placing controls and instruments at a comfortable level.


Dual level fluid sensors alert the user when the waste tank is 85% full for better procedure planning. The Transport III will then shut off the compressor before the tank exceeds capacity. The externally mounted waste container and optional amalgam separator means faster and cleaner waste disposal.

Transport III w/Scaler	AEU-525S	\$13,400
Transport III	AEU-525	\$11,700

Options		
Amalgam Separator Kit	AA-525AMAL	\$730
Micromotor w/ 6' cable	AE-240SC-40	\$1,550
Foot Switch with 10ft cable	AE-7PM-10	\$220
Replacement Parts		
Micromotor w/5' cable	AE-240SC-40	\$1,550
Waste Tank Assembly	330603	\$660
Foot Switch	AE-7PM	\$200
Power Cord 110V	840049	\$42
Power Cord 220V	840007	\$13
Water Reservoir (white)	730471	\$13
Replacement Scaler Tips For AEU-525S Only		
Slim Universal	ASC-10-PE38	\$130
Thin Subgingival	ASC-10-PE37	\$130
Power Universal	ASC-10-PE39	\$130



Handpieces sold separately
Recommended handpieces see pages 16-18.



Delivery module and user-interface panel mount externally, placing controls and instruments at a comfortable level. AEU-525 has water connection for external scalers. AEU-525S includes built in piezo ultrasonic scaler with LED illumination and three tips.

Product shown with high volume suction and saliva ejector suction options added



Lightweight Portable Delivery Systems

ADU-04ST

- Compact design and case under 15lbs (6.8kg)
- Customizable handpiece (high/slow speed) and vacuum options
- 3-way air/water syringe
- 25 ft. air line to connect to any air source
- Instrument tray
- Padded carrying case
- Carrying case: 25" x 11.5" x 11.5" (63.5 x 29 x 29cm)
- Requires air compressor for operation
- Quick disconnect options, .25 QD-Standard or 3/8 Schrader (optional)

Restorative System w/ 2 HPC Connections	ADU-04ST	\$1,990
---	----------	---------

Optional Upgrades		
High Volume Suction	AA-21ST	\$760
Saliva Ejector	AA-25ST	\$700
5-hole Fiber Optics	AA-19LED-04	\$720
6-pin Fiber Optics	AA-19A-06	\$720
Fiber Optics for 6-pin LED handpieces	AA-19EGEN-04T6*	\$260

*Will not function with incandescent bulbs. Consult your Aseptico sales rep for details.

Airpac III Portable Compressor

AA-100-110

- 1.0 horsepower, 100 psi oil-less motor
- Air delivery 3.1 CFM at 40 PSI
- Built-in wheels and pull handle
- Available in 110V
- Not available in Canada
- Weight: 46 lbs. (20.87 kg)
- Size: 17.93" W x 22.06" L x 10.43" H (45.5 cm W x 56.0 cm L x 26.5 cm H)

The Airpac III Portable Compressor is built by Aseptico to stand-up to heavy use and travel, the compressor is secured in its own super duty carrying case with built-in wheels and pull handle. The Airpac III includes an air storage tank and automatic pressure switch control.

110V Compressor	AA-100-110	\$2,250
-----------------	------------	---------



ASEPTICO PORTABLE SYSTEMS ARE MADE IN THE USA. CALL TODAY! 425.408.9347



AMC-25 Mobile Dental Cart

The **AMC-25 Mobile Cart** is a superior solution for providing general dental care in nursing homes, hospitals, or anywhere a mobile dental unit is needed. Available optional accessories allow for customizing this dental cart to meet your needs, including the ability to add electric handpieces, handpiece illumination, and much more!

- No plumbing required!
- 2 automatic high speed air handpiece connections (4-hole; non-fiber optic)
- 1 automatic low speed handpiece connection (4-hole; non-fiber optic)
- High and low volume suction with automatic overflow protection and solids trap
- 3-way air/water syringe
- Quick connect for ultrasonic scaler
- 2 hospital-grade outlets for accessories
- Wet/dry foot control
- Air Auxiliary Quick Disconnect
- Dual voltage 110V or 220V
- Operating pressure 80 psi
- High volume vacuum 7.6 SCFM @ 1.3 in Hg, 5.6 SCFM @ 4 in Hg
- Low volume vacuum 2.2 SCFM @ 1.5 in H
- Noise Level 65 dBA or less @ 3'4" (1 meter)
- Weight 167 lb (75.8 kg)
- Dimensions 23.5" W x 30" L x 36.5" H (56.7 x 76.2 x 92.7 cm)
- Optional: integrated electric motor, fiber-optic illumination for handpieces, ultrasonic scaler, curing light, amalgam separator, external vacuum/air/water hospital connections, auxiliary arm with accessory tray



AMC-25 Mobile Dental Cart		AMC-25	\$15,600
Options			
Integrated Electric Motor	330674		\$1,890
Valo Curing Light	ALU-25CLK-VALO		\$1,795
2 Fiber Optic HP Connections			
5-hole (ISO-B)	AA-19LED-25CK5		\$1,400
6-pin (ISO-C)	AA-19A-25CK		\$1,350
Piezo Scaler	ASC-10-25CKL		\$1,850
Amalgam Separator	AA-20AMAL-CK		\$880
Auxiliary Arm with Accessory Tray	AA-20A		\$520
Hospital Kit	AA-25HK		\$1,670
(external air, water and vac connections)			
Replacement Parts			
Micromotor	AE-240SC-40		\$1,550
Wet/Dry Foot Control	AA-43W		\$110
Air/Water Syringe	TA-90D		\$85
1L Water Reservoir – Clear	730631-01		\$70
1L Water Bottle – Opaque (each)	730471		\$13
Strainer	730745		\$34
Amalgam Separator Filter	730595		\$160
Foot Control Caddy	AA-25CDY		\$45

Unparalleled Features & Customization

The **AMC-25 Mobile Cart** is not only the most advanced dental cart on the market, but it also offers a range of customizability to suit your specific needs.



Industry Leading VALO Curing Light
Easy connect/disconnect with its own individual holder.

Scaler Connection
Optional ultrasonic scaler connection can be hardwired in.

Multiple Fiber Optic Connections
Optional connections for your fiber optic handpieces.

Integrated Electric Motor
Read more about functionality on the AEU-5000 (pg. 9).

Completely Mobile for Easy Transport
Easy to wheel into trucks or vans, down hallways, in elevators, ER's and operating rooms.

Wet/Dry Foot Pedal
Stays in place and can be activated from any angle.

Auxiliary Arm with Accessory Tray
360 degree maneuverability to be able to operate from any angle.

Amalgam Separator
Filters out particulate matter from waste water.

Foot Control Caddy
Easily store the foot control and other accessories on your AMC-15/25

1L Water Reservoir & Bottle
Large capacity for minimal water recycling between patients.

Great Backup Unit for Offices & Hospitals

Impact Resistant
Durable, wheeled case with comfortable handles.

Hospitable Grade Outlets
For any additional accessories needed.



ASEPTICO MOBILE DENTAL SYSTEMS ARE MADE IN THE USA. CALL TODAY! 425.408.9347

AMC-15 Mobile Dental Cart

Designed for Easy and Affordable Mobility, the AMC-15 Mobile Dental Cart is equipped for general dentistry with two automatic high speed air connections, one low speed handpiece connection, internal vacuum and air, water supply, large waste tank with automatic overflow protection, and amalgam separator.

- No plumbing required!
- 2 automatic high speed air handpiece connections (4-hole; non-fiber optic)
- 1 automatic low speed handpiece connection (4-hole; non-fiber optic)
- High & low volume suction with automatic overflow protection and solids trap
- 3-way air/water syringe
- Quick connect for ultrasonic scaler
- Amalgam Separator
- Wet/dry foot control
- Air Auxiliary Quick Disconnect
- Available in 110V or 220V
- Operating pressure 80 psi
- High volume vacuum 6.0 SCFM @ 0 in Hg, 3.6 SCFM @ 4 in Hg
- Low volume vacuum 1.3 SCFM @ 1.5 in Hg
- Noise Level 63 dBA @ 3'4" (1 meter)
- Weight 139 lb (63 kg)
- Dimensions 23.5"W x 30" L x 36.5" H (56.7 x 76.2 x 92.7 cm)
- Optional: auxiliary arm with accessory tray

AMC-15 Mobile Dental Cart **AMC-15** **\$12,500**

Options

Auxiliary Arm with Accessory Tray	AA-20A	\$520
Single Highspeed 6-Pin	AA-19EGEN-04T6	\$260
Fiber Optic Line (Add Up To 2 Per Cart)		

Recommended Handpieces

Air-Driven *see pages 16-17*

Replacement Parts

3-Way Air/Water Syringe with Tip	TA-90D	\$85
Amalgam Separator Filter	730595	\$160
Clear 1-liter Water Bottles w/Level Indicators	730631-01	\$70
Two 1-liter Water Bottles	730471	\$13
Wet/Dry Foot Control	AA-43W	\$110
Power Cord 115V	840049	\$42
Power Cord 230V	840007	\$13
Foot Control Caddy	AA-25CDY	\$45



"We have been using our AMC-15 since November 2023 approximately once per month for dental care in people with severe developmental disabilities. It is the perfect unit for our needs and has functioned flawlessly.

Gregory A. Costello, MD
Chief Medical Officer
Hope Christian Services



Compricart II Delivery System

ADU-20B

- One high speed and one low speed handpiece connection
- High volume evacuator and saliva ejector
- Large stainless steel instrument tray (13.5" x 9.75")
- All day water system with two liter self-contained water bottle

Roll the Compricart II Dental Cart from room to room with ease! Great for hospitals and expanding dental offices, the simple controls of the Compricart II make operation stress-free. A compact design, ample work surface, stainless steel instrument tray, integrated water supply, and venturi suction (high volume and saliva ejector) round out the features of this versatile platform. Just connect to an air source – no central vacuum required! Including the optional AWS-1C allows an external scaler to use the cart's water supply.



Compricart II	ADU-20B	\$3,600
Water Connection for External Scalers	AWS-1C	\$85

Options

Fiber Optic Upgrade Kit 5-hole	AA-19LED-04	\$720
Fiber Optic Upgrade Kit 6-pin	AA-19A-06	\$720

Replacement Parts

Complete HVE Assembly	330461	\$520
Disc Foot Control w/ Water Toggle	AA-43W	\$110
Air/Water Syringe	TA-90D	\$85
2L Water Bottle	730665	\$16



AEU-5000 Add-on Option Available
See details on page 9

FRIENDLY AND PROFESSIONAL SERVICE. CALL TODAY! 425.408.9347



NEW Sileather Fabric

Discover Sileather™ — eco-friendly and sustainable silicone fabric designed for durability and easy care. Naturally stain-resistant, UV-stable, and tough against wear, Sileather keeps our stools and chairs reliably clean and comfortable in demanding settings.



Chairs and Stools Available in Gray and Blue



All seating includes a padded carrying case



AseptiStool

ADC-08

- Set-up and store in seconds
- Includes padded case – no loose parts to misplace!
- 300 lb weight capacity
- Adjustable seat height: 19" - 22" (48.3-55.9 cm)
- 15 lbs (6.75 kg)
- Padded carrying case: 18" x 22" x 8" (48 x 53 x 23 cm)

Durable, adjustable and foldable, the AseptiStool Portable Dental Stool is the perfect seating solution for any clinician. Its one-piece design sets up or stores in seconds, and it swivels on four casters for easy mobility during treatment. It is height adjustable, lightweight, and folds up for simple storage in the included padded travel case.

AseptiStool	ADC-08	\$970
Replacement Padded Carrying Case	410093	\$180



Hydraulic AseptiStool

ADC-10

- Easy set-up and storage
- Includes padded case
- 300 lb weight capacity
- Hydraulic adjustable height: 18" - 25.5" (45.7 - 64.8 cm)
- 19 lbs (8.6 kg)
- Padded carrying case: 26" x 26" x 11" (56 x 56 x 28 cm)

With hydraulic adjustable height, the Hydraulic AseptiStool gives the mobile dentist the comfort and flexibility of doing dentistry in the office. It is easy to set up, tear down and transport in its included padded case. Even with its sturdy construction and high weight capacity, the stool is lightweight and carrying it into any mobile dental scenario is simple.

Hydraulic AseptiStool	ADC-10	\$1,300
Replacement Padded Carrying Case	410213	\$240



AseptiChair

ADC-01

- Sturdy, high-strength steel construction
- Lightweight, only 30 lbs (13.5 kg) in case!
- 300 lb weight capacity
- Backrest adjusts from fully upright to supine
- Price includes padded case, arm slings, and headrest cushion
- Adjustable seat height: 11" - 22" (11 - 55.9 cm)
- Padded carrying case: 23" x 34" x 8" (58 x 86 x 20.3 cm)

The AseptiChair's sturdy rolled steel construction is lightweight yet strong, making it perfect for portable dentistry. It is adjustable to five positions, with seat height ranging from 11 to 22 inches. The toe board folds down, allowing this portable dental chair to be used for examination.

AseptiChair	ADC-01	\$2,200
Pediatric AseptiChair (Red)	ADC-01P-RED	\$2,200

Replacement Parts		
Replacement Padded Carrying Case	410193	\$280
Replacement Leg Adjustment Pin	460445-08	\$17
Endcap Heavy Duty 1"	850026	\$3
Armrest Support Pin	460444-08	\$9

"Having quality products like Aseptico patient/operator chairs gives our providers, patients, and facilities confidence in the care that a mobile alternative practice can provide"

Sara Luetke, RDH
Sound Dental Care



Hydraulic AseptiChair

ADC-02

- Hands-free height and back adjustment
- Adjustable with patient in chair
- 300 lb weight capacity
- No external electricity or air required
- Handle and wheels for easy mobility
- 69 lbs (31.3 kg)
- Hydraulic adjustable seat height: 11" - 33" (27.9 - 83.8 cm)
- Folded dimensions: 33" x 22" x 18" (83.8 x 55.9 x 45.7 cm)

The Hydraulic AseptiChair provides hands-free operation and is the only portable chair available with adjustability functions comparable to a standard dental chair. All chair lift functions are foot actuated hydraulic for hands-free operation; no external power is required. The AseptiChair provides comfortable patient seating for patients up to 300 lbs. Easy move wheels.

Hydraulic AseptiChair	ADC-02	\$7,250
Hydraulic Assembly Only	330629	\$2,000

Replacement Parts		
Replacement Dust Cover	410210	\$420



ALL ASEPTICO SEATING IS MADE IN THE USA. CALL TODAY! 425.408.9347

NOMAD™ PRO 2 Defining Portable Dental X-ray Technology

- ARU-07KP2
- Best-in-class portable digital X-ray system
 - Lightweight 5.5lbs (2.5kg)
 - Simple touch controls on color LCD screen
 - Removable and rechargeable battery packs
 - One battery charge allows for more than 100 exposures
 - Heavy Duty Protective Transport Case Included

The NOMAD Pro 2, Portable Dental X-ray System captures clear radiographic images with ease – anytime, anywhere. Its lightweight, ergonomic, handheld design makes it ideal for mobile dental units, or practices with multiple operatories. Best-in-class features include advanced shielding to ensure operator and patient safety while delivering precise imaging results.



Offering cordless, handheld operation, the NOMAD Pro 2 portable x-ray system is powered by rechargeable battery packs and comes contained in an easily transportable, rugged case with pull-out handle and wheels. Long-lasting batteries deliver exceptional performance, allowing you to complete hundreds of exposures on a single charge. Intuitive user interface and quick setup save valuable time, enabling you to focus on what matters most—providing patient care.

Nomad Pro 2 Portable X-Ray System with Carrying Case	ARU-07KP2	Call For Best Price
Includes Premium Upgrade Kit		
Replacement Parts		
Premium Upgrade Kit		
Includes Two Replacement Battery Handsets		
One Replacement Charging Cradle		
Power Cord		
	ARU-10134502	\$953
Basic Upgrade Kit		
Includes Replacement Battery Handset		
One Replacement Charging Cradle		
Power Cord		
	ARU-10135062	\$589
Replacement Battery Handset	ARU-10131887	\$519
Replacement Charging Cradle		

Not Available In Canada



Portable LED Exam Light

- ALU-40LED
- Adjustable LED with seven beam diameters
 - Durable construction
 - Includes padded carrying case
 - Optional Amber UV filter for working with composites
 - 12 lbs (5.4 kg)
 - Padded carrying case: 45" x 9" x 5.5" (114.5 x 23 x 14 cm)
- The Portable LED Exam Light** provides high-intensity LED light with seven beam diameter settings. It features a compact lamp head with adjustable intensity, flexible gooseneck, telescoping tripod stand and padded carrying case. When a standard electric outlet is not available, an optional auxillary DC input can be purchased that allows 12V and 24V battery operation in the field.

Portable LED Exam Light	ALU-40LED	\$1,600
Optional Amber UV Filter	461974	\$4
Replacement Padded Carrying Case	410084	\$170
DC Cable with Battery Clips	AA-40DC-BC	\$150
DC Cable with Car Lighter Adapter	AA-40DC-LP	\$170


Portable Tray Stand

ATC-03CF

The Portable Tray Stand is a sturdy and necessary addition to any mobile operatory. Standard size stainless instrument tray removes easily for autoclaving.

It weighs less than 9 pounds and folds into included case when work is done. Made in USA.

Portable Tray Stand	ATC-03CF	\$990
---------------------	----------	-------



CALL TODAY! 425.408.9347
FOR FRIENDLY AND PROFESSIONAL CUSTOMER SERVICE.

Heavy Duty
Brushless Lab Motor
Torque Pro
Delivers Power

AEU-16SS-120

- For heavy acrylic use
- Brushless 50k rpm micromotor
- Consistent torque even at low speeds
- Includes Variable Speed Foot Control



The Heavy Duty Brushless Lab Motor delivers power, performance, durability, and portability. The brushless 50k rpm micromotor can produce up to 8.3 Ncm of torque and features easy to read digital display, preset memory storage, and hand/foot change button.

Heavy Duty Dental Lab Motor **AEU-16SS-120** \$1,150

Replacement Parts:

Motor/Cable Assy	AE-H180	\$540
Variable Speed Footswitch	AE-SP	\$75

Light Duty Lab Motor
Powerful & Lightweight

AEU-10SS

- For light to medium acrylic use and low speed proph
- Versatile E-type straight 1:1 handpiece included
- 20k rpm micromotor
- Available in 110V or 220V
- Includes AE-7P On/Off Foot Switch

The Light Duty Lab Motor is one of the finest motor and handpiece combinations available to the professional. It is ideal for chairside adjustments, prophylaxis, or any light to medium use.

Light Duty Dental Lab Motor **AEU-10SS** \$1,150

Light Duty Dental Lab Motor **AEU-10SS-220** \$1,150

Options

Optional Variable Speed Foot Control	AE-14A	\$360
--------------------------------------	---------------	-------

Replacement Parts:

1:1 Handpiece	AHP-101	\$380
Micromotor	AM-103B	\$220
On/Off Foot Switch	AE-7P	\$190

The Battery Powered Portable Lab Motor is a smart solution for lightweight lab benchwork and performing prophylaxis in the field! This small portable unit delivers 22k rpm and can perform 8 hours on a single battery charge.

Portable Lab Motor with Handpiece **AEU-03M1B** \$800

Portable Lab Motor **AEU-03M1B-NHP** \$570

Replacement Parts

1:1 Handpiece	AHP-101	\$380
Micromotor	AE-5M1-25	\$180
Foot Switch	AE-7M1	\$23

Battery Powered
Portable Lab Motor
Ideal for Mobile Hygiene

AEU-03M1B

- For light duty use and low speed proph
- Straight 1:1 e-type handpiece included
- Dual voltage 110V-240V



Inludes AHP-101 Handpiece

Shop online at aseptico.com

Portable Sterilizer

ASU-12-110 and ASU-12-230

- Compact and portable
- 18-minute sterilization cycle (24-minute total cycle time)
- Only 6 minutes to reach sterilization temp
- Temperature: 250°F (121°C)
- 11.5 lbs (5.2 kg)
- 9 liter capacity

The Portable Dental Autoclave is designed to sterilize solid, unwrapped instruments in a 121°C saturated steam process during 18-minute sterilization cycle. Use distilled or deionized water only.

110V Sterilizer **ASU-12-110** \$1,399

230V Sterilizer **ASU-12-230** \$1,499



Portable Saliva Ejector

AA-20

- Great for mobile dentistry or as a backup for sedation dentistry
- Powered by rechargeable battery pack, wall outlet or 12 VDC vehicle source
- Up to 60 minutes of continuous use on a single charge
- 800cc waste canister with float shut-off
- 6.3 lbs (2.9 kg)
- 9" x 7" x 8" (22.8 x 17.8 x 20 cm)

Great for mobile dentistry or as a backup for sedation dentistry, the Portable Saliva Ejector can be powered by a wall outlet, 12V DC source, or rechargeable battery pack. It features an adjustable 80 - 550 mm Hg vacuum range and 800cc waste canister with float shut-off. Includes padded carrying case.

Portable Saliva Ejector **AA-20** \$620

Replacement 800 ml Canister **AA-20CAN** \$38

Bacteria Filters **AA-20BF-2** \$12

Vacuum Pump Tube Assembly **330668** \$49



In addition to providing easy portable suction, the Portable Saliva Ejector fulfills the legal requirement for a back-up suction system in the event of power or suction failure.



CALL TODAY! 425.408.9347
FOR FRIENDLY AND PROFESSIONAL CUSTOMER SERVICE.

Customer Service

425.408.9988

Monday-Friday
6:00 a.m.- 3:30 p.m. PST
Customer service representatives
are standing by to help.
repairs@aseptico.com

Quality System

Aseptico's Quality Management System is certified to ISO 13485 as well as the Medical Devices Single Audit Program Recognized By The U.S. and Canada.

Return and Repair Policies

- Products that have been in clinical use cannot be returned for credit.
- Service and repair charges are billed on credit card. Customer is expected to accept charges up to \$750 without advance notification. Minimum service/repair charge is \$160.
- All returns and repairs require prior authorization. Contact Aseptico Customer Service for a Return Material Authorization number (RMA#).
- The RMA# must appear clearly on the exterior of the package being returned or it may be rejected and returned to the sender.
- Packages received without prior authorization are not the responsibility of Aseptico. Aseptico reserves the right to dispose of such packages and their contents 30 days after receipt with no obligation to the sender.
- The sender is responsible for proper packaging and for all costs for return shipment to Aseptico.
- New or unused products may be returned for full credit (less shipping charges) under a Return Material Authorization number within 30 days of date of original invoice.

Warranty

Aseptico warrants its products against defects in material or workmanship for a period of two (2) years, from date of original invoice. This warranty is not transferrable. Some handpieces are warranted for one year under the same conditions. Other handpieces and expendable components, such as air turbines and light bulbs, are covered by shorter warranty periods, or have no warranty. Aseptico's sole obligation under product warranty is (at its sole option and discretion) to repair or replace any defective component or product in part or whole. Aseptico shall be the sole arbiter of such action. In the event of alleged defect under warranty, the purchaser is to notify Aseptico's Customer Service Department promptly. Customer Service will provide instructions, usually directing that the product be returned for service. Shipment to Aseptico and the cost thereof is always the responsibility of the purchaser.

Accidental misuse, inappropriate installation, or failure to perform directed maintenance voids the warranty.

Aseptico does not assume, under this warranty, any risks or liabilities arising from the clinical use of its products, whether or not such use involves coincidental utilization of products manufactured by others. Deliberately defacing, modifying, or removing the serial number voids the warranty.

Sales

425.408.9347

Monday-Friday
6:00 a.m.- 4:30 p.m. PST
Sales representatives
are ready to assist you.
orders@aseptico.com

Ordering Policies

- Aseptico accepts Visa, MasterCard, Discover, and American Express.
- Open account is available only if credit is approved in advance.
- Freight costs are added to all invoices.
- All orders are shipped FOB origin.
- Delinquent accounts are subject to 1.5% monthly interest, attorney's fees, collection and court costs when needed to satisfy debt.
- Prices subject to change without notice.
- Please visit: aseptico.com/warranty-policies for additional warranty and repair policy information.

Sales Tax

Aseptico uses a tax service (Avalara) to calculate and apply state sales tax to invoices for all states that Aseptico is registered to collect/remit sales tax. We reserve the right to register in additional states at any time and begin collecting/remitting sales tax.

The amount of tax charged on your order will depend upon many factors including type of item purchased and destination of the shipment. Factors can change between the time you place an order and the time of invoicing, which could affect the calculation of sales taxes. The amount appearing on your order as sales tax may differ from the sales taxes ultimately charged.

Customer is responsible for any applicable state sales or use taxes not charged by Aseptico.

Sales Tax Exemptions

Unless Aseptico is in possession of a validated exemption certificate prior to shipment, the customer must pay the entire invoice, including sales tax. Certificates can be emailed to AR@aseptico.com. You can upload your Tax Exempt certificate using our Avalara Cert Express Tool: aseptico.com/tax-exempt. Please inform your sales rep when you do so, and they will get it applied to your account for future orders placed through your sales rep.

International Customers

Most Aseptico products are CE-marked and Health Canada licensed. However, customers should verify the status of any specific item(s) with their sales representative or local distributor prior to finalizing the purchase. Please contact your local Aseptico distributor for pricing and availability.

Prices shown in the catalog are in US dollars. If you are unsure of the Aseptico distributor in your country, for specific ordering and payment policies, or for general product inquiries, please contact our International Sales Department at int@aseptico.com.

Aseptico Product Recognition

Aseptico® Products Have Been Recognized Multiple Times by leading dental industry publications for their innovation and preference for use by dental professionals. *Dentistry Today*, the leading US clinical magazine for dentists, recently recognized multiple Aseptico products including:



GO Ultra Portable AEU-350S

- July/Aug 2024 – Dentistry Today Top 50 Technology Products



AMC-15 Mobile Dental Cart

- Jan/Feb 2024 – Dentistry Today Reader's Choice Top Innovative Products



Aseptico® 7000 Series All-in-One Motor

- March 2024 – Dentistry Today Reader's Choice Top 25 Implant Product
- April 2024 – Dentistry Today Reader's Choice Top 25 Endodontic Product
- October 2024 – Dentistry Today Editor's Choice Awards
- Top 25 Implant Product and Top 25 Endodontic Product for multiple years running!



1070 Series Implant Motor

- July/Aug 2024 – Dentistry Today Top 100 Products

AEU-1070-70V & AHP-85-SS

1070 Motor With FREE Handpiece Package

Only \$3,600

Receive a FREE AHP-85-SS implant handpiece – a value of \$670 when you purchase the New 1070 Series Implant Motor

Call 425-408-9347 To Learn More

NEW 20:1 Economy
Dental Implant Handpiece
AHP-85-SS

1070 Series
Implant Motor
AEU-1070-70V



RFA-55250-MULTI

PENGUIN II

Call for Best Price Available

- Accurate ISQ measurement for greater confidence determining when to load the implant
- Cordless lightweight and simple to use
- Half the price of competing ISQ devices

Call 425-408-9347 To Learn More

